

STERRAD® 100NX® SYSTEM CYCLE SELECTION



**The Next
Generation
Sterilization
System.**

All instruments must be cleaned and thoroughly dried prior to being packaged for sterilization in the STERRAD® 100NX® Sterilizer.

Approved packaging includes ASP SteriTite® Containers, wrapped or pouched APTIMAX® Trays or Tyvek® Pouches and Rolls with STERRAD® Chemical Indicator Strip. For rigid containers, verify that they are cleared for use in each cycle.

ADVANCED STERILIZATION PRODUCTS

Division of Medos International Sàrl
a **Johnson & Johnson** company

*Please refer to the STERRAD® 100NX® Sterilization System User's Guide for important information prior to use. Consult with the medical device manufacturer of instruments or endoscopes you intend to sterilize with any questions regarding materials compatibility.

All Rights Reserved.
Tyvek and Teflon are registered trademarks of E.I. du Pont de Nemours and Company.
da Vinci is a registered trademark of Intuitive Surgical, Inc.

STANDARD CYCLE: 47 MINUTES

The **STANDARD** Cycle should be selected for instruments that meet the criteria below:

- General medical instruments (metal and nonmetal, including hinged devices)
- Instruments with single-channel stainless steel lumens
 - Internal diameter 0.7 mm or larger and length 500 mm or shorter
- Polyethylene and Teflon lumen tubing with an internal diameter of 1 mm or larger and length of 1000 mm or shorter.
 - These items must be sterilized without any additional load items
 - Limit of 20 pieces of tubing per cycle
- **Maximum of 40 single channel stainless steel lumens**

Instrument sets that can be sterilized in the **STANDARD** Cycle include but are not limited to*:

- Arthroscopy and laparoscopic instrument sets
- Eye instruments
- Cystoscopy instruments
- Rigid or semi-rigid ureteroscopes
- Cameras and light cords
- Rechargeable batteries
- Doppler cords and defibrillator paddles
- Orthopedic drills and saws
- Ultrasound probes/transducers



STANDARD

Load Preparation: Two shelves (total weight: 21.4 lbs or 9.7 kg)

FLEX CYCLE: 42 MINUTES

The **FLEX** Cycle should be selected for flexible endoscopes that meet the criteria below:

- Single channel
- Internal lumen diameter of 1 mm or larger and length of 850 mm or shorter
- A maximum of two single-channel flexible endoscopes may be processed at one time
- No additional items may be processed in the cycle with the flexible endoscopes
- Place venting caps on flexible scopes according to the manufacturer's instructions

Flexible endoscopes that can be sterilized in the **FLEX** Cycle include but are not limited to*:

- Bronchoscopes
- Cystoscopes
- Hysteroscopes
- Choledochoscopes
- Flexible ureteroscopes
- Thorascopes
- Intubation fiberscopes



FLEX

Load Preparation: A maximum of two single-channel flexible endoscopes

EXPRESS CYCLE: 24 MINUTES

The **EXPRESS** Cycle should be selected for instruments that meet the criteria below:

- General metal medical devices requiring surface sterilization, or sterilization of mated stainless steel and titanium surfaces
- Rigid or semi-rigid endoscopes without lumens
- Rechargeable batteries

Materials that should not be processed in the **EXPRESS** Cycle, even though they can be processed in the **STANDARD** and **FLEX** Cycles, are:

- Items made of nylon, polyurethane, or Kraton®
- Items with mated Delrin®, mated Ultem®, mated Radel®, or mated aluminum surfaces

Instruments that can be sterilized in the **EXPRESS** Cycle include but are not limited to*:

- da Vinci® endoscopes
- Rigid or semi-rigid endoscopes without lumens
- General surgery metal devices without lumens
- Rechargeable batteries
- Eye Instruments without lumens



EXPRESS

Load Preparation: Bottom shelf only (10.7 lbs or 4.85 kg)