



NU-CIDEX®

PERACETIC ACID SOLUTION

HIGH-LEVEL DISINFECTION

NU-CIDEX® Peracetic Acid Solution has been designed by ASP to high-level disinfect flexible endoscopes and other heat-sensitive devices that cannot be terminally sterilized.

EXPERIENCE THE BENEFITS

NU-CIDEX® Peracetic Acid Solution contains 0.35% peracetic acid and is cleared as a high-level disinfectant.

- **Effective with sporicidal action**
- **Fast-acting**
- **Safer for medical personnel and environment**





NU-CIDEX®

PERACETIC ACID SOLUTION

FEATURES & BENEFITS

- **EFFECTIVE WITH SPORICIDAL ACTION** – achieves high-level disinfection and sporicidal activity in only 5 minutes at room temperature.
- **FAST-ACTING** – fast turnaround of reprocessed instruments.
- **SAFE TO USE FOR HEALTHCARE PROFESSIONALS** – low vapor pressure of in-use solution.
- **SAFER FOR ENVIRONMENT** – naturally degrades into biodegradable components.

MATERIALS COMPATIBILITY

NU-CIDEX® Peracetic Acid Solution can be used safely to disinfect devices made from metals, plastics and elastomers. The solution contains corrosion inhibitors for an increased tolerance on brass and copper instruments that can be soaked for the recommended contact time.

MICROBICIDAL ACTIVITY

NU-CIDEX® Peracetic Acid Solution provides a wide spectrum activity against bacteria, mycobacteria, viruses, fungi and bacterial spores (sporicidal at 6 log reduction).

TECHNICAL INFORMATION

MEDICAL DEVICE CLASSIFICATION	Class IIb according to MDD 93/42/EEC
IN-USE CONCENTRATION	0.35% peracetic acid
SOAK TIME	High-level Disinfection: 5 minutes at room temperature Sporicidal activity (6 log reduction): 5 minutes at room temperature
USE LIFE	24 hours
SHELF LIFE	35 weeks unopened
OPEN BOTTLE SHELF LIFE	24 hours (when open and mixed)
STORAGE	Between 5°C-28°C, upright
DISPOSAL	Drain or as per hospital policy. Flush thoroughly with water

ORDERING INFORMATION

REORDER NO.	DESCRIPTION	CASE CONTENTS
NCX010	NU-CIDEX® Peracetic Acid Solution	2 x 5 liters



NU-CIDEX® Peracetic Acid Solution may also be used for instrument reprocessing in automated equipment, in accordance with equipment manufacturers' recommendations.

Please read and follow the Instructions for Use prior to using NU-CIDEX® Peracetic Acid Solution for detailed information, including contraindications, warnings and proper directions for use.

FOR MORE INFORMATION, CONTACT YOUR LOCAL ADVANCED STERILIZATION PRODUCTS SALES REPRESENTATIVE

www.aspjj.com/emea

<p>UK: Advanced Sterilization Products A Division of J&J Medical Ltd. Pinewood Campus Nine Mile Ride Wokingham Berkshire RG40, 3EW, England T: +44 1344 871 081 F: +44 1344 871 171</p>	<p>Ireland: Advanced Sterilization Products A Division of J&J Medical Ireland Airton Road, Tallaght, Dublin 24, Ireland T: +353 1 466 5200 F: +353 1 466 5340</p>	<p>Egypt: Johnson & Johnson Medical Egypt Florida Mall 5th Floor 1229 Square El Sheikh Ali Gad El Hak St. Heliopolis Cairo, Egypt T: +202 2268 5026 F: +202 2268 4674</p>	<p>Middle East: Advanced Sterilization Products A Division of J&J Medical Middle East Dubai Healthcare City, J&J building, 3rd floor PO Box 505080, Dubai, United Arab Emirates T: +971 4 429 7235 F: +971 4 429 7250</p>	<p>South Africa: Advanced Sterilization Products A Division of J&J Medical Pty. Ltd. SA PO Box 273, Midrand Halfway House 1685 South Africa T: +27 11 265 1120 F: +27 11 265 1189</p>
--	--	--	--	--

ADVANCED STERILIZATION PRODUCTS

Division of Medos International Sàrl

a **Johnson & Johnson** company