



Endoscope
Process Wheel

Pure Technology



ASP
CIDEX® OPA•C

ASP
ADAPTACIDE™ PAA-C


As an integral part of Advanced Sterilization Products' Total Solution for Endoscope Reprocessing, ASP offers a choice of disinfectants depending on user preference.

For high-level disinfection, CIDEX® OPA-C Solution is the recommended choice, however when sporicidal activity is sought, ADAPTACIDE™ PAA-C Solution is the preferred option.



Protecting Lives Against Infection™

Think solution, not product



ASP® CIDEX® OPA-C

CIDEX® OPA Concentrate Solution provides a fast and effective solution for high-level disinfection of flexible endoscopes. Based on the clinically proven CIDEX® OPA Solution, CIDEX® OPA-C Solution is a concentrated high-level disinfectant designed for use in the ADAPTASCOPE™ single-shot Washer-Disinfector (WD).

Concentration

5.75%, diluted to an in-use concentration of 0.0575%

Shelf Life

15 months

Minimum effective temperature

51°C

Medical Device Classification

Class IIa according to MDD 93/42/EEC

Disposal

Drain or as per hospital policy

Storage information

Store at controlled room temperature 15°C - 30°C

Product code

20396

Sales Unit

4 x 5 litres

CIDEX® OPA-C Solution is designed for use in the ADAPTASCOPE™ System. For use in other single-shot automated endoscope reprocessors, please contact your local ASP representative.



ASP® ADAPTACIDE™ PAA-C

ADAPTACIDE™ PAA-C Solution provides a fast and effective sporicidal solution for disinfection of flexible endoscopes, designed for use in the ADAPTASCOPE™ single-shot Washer-Disinfector (WD).

Concentration

15%, diluted to an in-use concentration of 0.15%

Shelf Life

12 months

Minimum effective temperature

25°C

Medical Device Classification

Class IIa according to MDD 93/42/EEC

Disposal

Drain or as per hospital policy

Storage information

Store at controlled room temperature 0°C - 25°C

Product code

20514

Sales Unit

2 x 5 litres

ADAPTACIDE™ PAA-C Solution is designed for use in the ASP ADAPTASCOPE™ System. For use in other single-shot automated endoscope reprocessors, please contact your local ASP representative.



Contact for UK:

ASP c/o J&J Medical Ltd.
Pinewood Campus, Nine Mile Ride
Wokingham
Berkshire RG40 3EW, England
T : +44 1 344 871 081
F : +44 1 344 871 171

Contact for Ireland:

ASP c/o J&J Medical Ireland
Airton Rd, Tallaght
Dublin 24, Ireland
T : +353 1 466 5200
F : +353 1 466 5340

Contact for Egypt:

Johnson & Johnson Medical Egypt
Florida Mall 5th Floor
1229 Square El Sheikh Ali Gad
El Hak St. Heliopolis
Cairo, Egypt
T : +202 2268 5026
F : +202 2268 4674

Contact for Middle East:

ASP c/o J&J Medical Middle East
Dubai Healthcare City,
J&J building, 3rd floor
PO Box 505080, Dubai,
United Arab Emirates
T : +97 1 4 429 7235
F : +97 1 4 429 7250

Contact for South Africa:

ASP c/o J&J Middle Pty Ltd SA
PO Box 273, Midrand
Halfway House 1685
Republic of South Africa
T : +271 1 265 1120
F : +971 1 265 1189

ADVANCED STERILIZATION PRODUCTS

Division of Medos International Sàrl

a **Johnson & Johnson** company