



STERRAD® SYSTEM STERILIZATION RECORD
STERRAD® 100NX® System

DAILY BIOLOGICAL INDICATOR (BI) TEST CONTROL SHEET

TODAY'S DATE:		BIOLOGICAL INDICATOR (BI) LOT NUMBER(S):	ST - FL - EX- Control -
TIME:		Incubator Temperature BETWEEN 55°C and 60°C	

24-HOUR BIOLOGICAL INDICATOR (BI) STANDARD	24-HOUR BIOLOGICAL INDICATOR (BI) FLEX	24-HOUR BIOLOGICAL INDICATOR (BI) EXPRESS
Date Read:	Date Read:	Date Read:
Time Read:	Time Read:	Time Read:
Read by:	Read by:	Read by:

BI Negative (yellow cap and purple media: Processed BI)	Results - +	BI Negative (yellow cap and purple media: Processed BI)	Results - +	BI Negative (yellow cap and purple media: Processed BI)	Results - +
CONTROL Positive (red cap and yellow media: Positive BI)	- +	CONTROL Positive (red cap and yellow media: Positive BI)	- +	CONTROL Positive (red cap and yellow media: Positive BI)	- +

Indicate type of Cycle: **STERRAD® 100 NX® System**
 STANDARD FLEX EXPRESS

CYCLE NUMBER	CYCLE LOAD CONTENTS	NO. OF ITEMS	CHEMICAL INDICATOR AND/OR TAPE		INITIALS
			PASS	FAIL	

Indicate type of Cycle: **STERRAD® 100 NX® System**
 STANDARD FLEX EXPRESS

CYCLE NUMBER	CYCLE LOAD CONTENTS	NO. OF ITEMS	CHEMICAL INDICATOR AND/OR TAPE		INITIALS
			PASS	FAIL	

This sheet is designed for use with STERRAD® CYCLESURE® 24 Biological Indicator. Please refer to the complete Instructions For Use for the STERRAD® CYCLESURE® 24 Biological Indicator.

ADVANCED STERILIZATION PRODUCTS®

Division of Ethicon, Inc.

a **Johnson & Johnson** company

AD-110402-01-US_A