



**STERRAD® SYSTEM STERILIZATION RECORD**

**STERRAD® NX® System**

**DAILY BIOLOGICAL INDICATOR (BI) TEST CONTROL SHEET**

<b>TODAY'S DATE:</b>		<b>BIOLOGICAL INDICATOR (BI) LOT NUMBER(S):</b>	<b>ST - ADV - CONTROL -</b>
<b>TIME:</b>		<b>Incubator Temperature BETWEEN 55°C and 60°C</b>	

<b>24-HOUR BIOLOGICAL INDICATOR (BI) STANDARD CYCLE</b>
<b>Date Read:</b>
<b>Time Read:</b>
<b>Read by:</b>

<b>24-HOUR BIOLOGICAL INDICATOR (BI) ADVANCED CYCLE</b>
<b>Date Read:</b>
<b>Time Read:</b>
<b>Read by:</b>

<b>BI Negative</b> (yellow cap and purple media: Processed BI)	<b>Results</b> - +
<b>CONTROL Positive</b> Positive (red cap and yellow media: Positive BI)	- +

<b>BI Negative</b> (yellow cap and purple media: Processed BI)	<b>Results</b> - +
<b>CONTROL Positive</b> Positive (red cap and yellow media: Positive BI)	- +

Indicate type of Cycle: **STERRAD® NX® System**  
STANDARD  ADVANCED

CYCLE NUMBER	CYCLE LOAD CONTENTS	NO. OF ITEMS	CHEMICAL INDICATOR AND/OR TAPE		INITIALS
			PASS	FAIL	

Indicate type of Cycle: **STERRAD® NX® System**  
STANDARD  ADVANCED

CYCLE NUMBER	CYCLE LOAD CONTENTS	NO. OF ITEMS	CHEMICAL INDICATOR AND/OR TAPE		INITIALS
			PASS	FAIL	

This sheet is designed for use with STERRAD® CYCLESURE® 24 Biological Indicator. Please refer to the complete Instructions for Use for the STERRAD® CYCLESURE® 24 Biological Indicator.

**ADVANCED STERILIZATION PRODUCTS®**

Division of Ethicon, Inc.

a *Johnson & Johnson* company

AD-110401-01-US\_A