



CIDEZYME[®] XTRA

MULTI ENZYMATIC DETERGENT

Fast. Effective. Thorough.

CIDEZYME[®] XTRA Multi Enzymatic Detergent is formulated for the cleaning cycle in automatic endoscope reprocessors or manually. Containing multiple enzymes, it provides rapid, more complete penetration and effective removal of bioburden.

EXPERIENCE THE BENEFITS

CIDEZYME[®] XTRA Multi Enzymatic Detergent is an efficient first step before high-level disinfection

- Advanced, low-foaming solution is biodegradable and rinses easily
- Contains a rust inhibitor to help protect valuable medical devices and surgical instruments
- Easily and naturally removes fats, lipids, proteins, starches and carbohydrates





CIDEZYME[®] XTRA

MULTI ENZYMATIC DETERGENT

FEATURES

MULTI ENZYME FORMULATION

Cellulase, lipase, protease, amylase, and proprietary enzymes with advanced proteolytic action

LOW FOAMING

CORROSION INHIBITORS

BIODEGRADABLE

BENEFITS

Increased protein enzyme activity to provide a faster, more thorough penetration into hard-to-reach places

Effectively removes bioburden

For use in manual reprocessing and automated systems

Is compatible with medical devices

Safe for the environment

Materials Compatibility

CIDEZYME[®] XTRA MULTI ENZYMATIC DETERGENT HAS BEEN TESTED WITH A WIDE RANGE OF MATERIALS COMMON TO FLEXIBLE ENDOSCOPES INCLUDING:

Brass, copper, stainless steel, a variety of polymers, rubber, corrugated tubing and plastic

System Use

VALIDATED FOR USE IN THE ENDOCLENS-NSX[™] AUTOMATIC ENDOSCOPE REPROCESSOR (AER) DESIGNED FOR EITHER MANUAL OR AUTOMATIC REPROCESSING

FOR AUTOMATED SYSTEMS 4 ml per 1 l of water

FOR MANUAL USE 8 ml per 1 l of water

SHELF-LIFE 18 months

STORAGE INFORMATION 15-30°C

ORDERING INFORMATION

REORDER NO.	DESCRIPTION	CASE CONTENTS
22591	3.8 l container	2 x 3.8 l

Please read and follow the CIDEZYME[®] XTRA Multi Enzymatic Detergent instructions for use prior to using this product for complete usage information, including contraindications and warnings.

ENDOCLENS-NSX[™] is not CE marked and not sold in EU Countries.



FOR MORE INFORMATION PLEASE VISIT WWW.ASPJJ.COM/EMEA OR CONTACT YOUR ASP LOCAL REPRESENTATIVE

Egypt:

Johnson & Johnson Medical Egypt
Florida Mall
5th Floor
1229 Square El Sheikh Ali Gad El
Hak St. Heliopolis
Cairo, Egypt

T: +202 2268 5026
F: +202 2268 4674

Middle East:

Advanced Sterilization Products
A Division of J&J Medical Middle East
Dubai Healthcare City, J&J building,
3rd floor
PO Box 505080, Dubai,
United Arab Emirates

T: +971 4 429 7235
F: +971 4 429 7250

South Africa:

Advanced Sterilization Products
A Division of J&J Medical Pty.
Ltd. SA
PO Box 273, Midrand
Halfway House
1685
South Africa

T: +27 11 265 1120
F: +27 11 265 1189

ADVANCED STERILIZATION PRODUCTS

Division of Medos International Sàrl

a *Johnson & Johnson* company