



ENDOCLENS-NSX™

AUTOMATIC ENDOSCOPE REPROCESSOR

Proven. Consistent. Tested.

Introducing the ENDOCLENS-NSX™ Automatic Endoscope Reprocessor (AER), an automated high-level disinfection system that helps safeguard your people, patients and medical instruments while freeing up time for better care.

EXPERIENCE THE BENEFITS

Designed to promote safety and improve efficiency of endoscope reprocessing

- Over 5,000 systems in use with strong customer satisfaction
- Automation provides reliable results to promote a higher standard of care
- Integrated system offers complete high-level disinfection with world-leading CIDEX® OPA Solution






ENDOCLENS-NSX™

TECHNICAL SPECIFICATIONS

SYSTEM USE

TOTAL CYCLE TIME	19 minutes (14 minutes without wash stage)	
CAPACITY/BASIN	Basins per system	1
	Scopes per basin	1
	Basin size	47.0 cm depth x 39.5 cm width
	Basin volume	15 l (maximum)
	Basin Chamber Materials	Fiber reinforced plastic
DIMENSIONS AND HEIGHT	Height (lid closed)	96.5 cm
	Width	53.0 cm
	Depth	79.6 cm (66.5 without filter)
	Weight	74 kg (dry, including water filter)
ELECTRICAL REQUIREMENTS	ENDOCLENS-NSX™ AER Single Phase; 120/220/230/240 V ~; 0.8kVA; 50/60 Hz	
SCOPE SETTING	Loading	Top loading
	Scope connection required	Yes
	Number of channels	4
	Applicable Scopes	OLYMPUS/FUJINON/PENTAX
WATER REQUIREMENTS	Pressure	15 psi - 100 psi
	Hardness	< 60 ppm CaCO ₃
	Incoming Temperature	Ambient 5°C ~ 40°C
	Consumption per Cycle	55 liters
	Fill Rate	12 liters per minute
GERMICIDAL AGENTS	Approved Agents CIDEX™ OPA Solution	
OPERATIONAL ENVIRONMENT	Temperature	Ambient 20°C ~ 40°C
	Relative Humidity	Ambient 10-90% non-condensing)
	Elevation	400 m below sea level – 3000 m above sea level
FLOOR	Flat surface with no slope	Within 1 degree of vertical
FILTER MEDIA	0.2 micron filter	Polyether Sulfone or
		Polypropylene Polysulfon

Front View



Side View



ORDERING INFORMATION

CODE	PRODUCT DESCRIPTION
27000	ENDOCLENS-NSX™ AER
20391	CIDEX® OPA Solution (2 x 5 l/case)
22591	CIDEZYME® XTRA Multi Enzymatic Detergent (2 x 3.8 l/case)
20515	ADAPTACLEAN® Alkaline Detergent (2 x 5 l/case)
104714-01	0.2 Micron Filter

ENDOCLENS-NSX™ is not CE marked and not sold in EU Countries

FOR MORE INFORMATION PLEASE VISIT WWW.ASPJJ.COM/EMEA OR CONTACT YOUR ASP LOCAL REPRESENTATIVE

Egypt:

Johnson & Johnson Medical Egypt
Florida Mall
5th Floor
1229 Square El Sheikh Ali Gad El
Hak St. Heliopolis
Cairo, Egypt

T: +202 2268 5026
F: +202 2268 4674

Middle East:

Advanced Sterilization Products
A Division of J&J Medical Middle East
Dubai Healthcare City, J&J building,
3rd floor
PO Box 505080, Dubai,
United Arab Emirates

T: +971 4 429 7235
F: +971 4 429 7250

South Africa:

Advanced Sterilization Products
A Division of J&J Medical Pty.
Ltd. SA
PO Box 273, Midrand
Halfway House
1685
South Africa

T: +27 11 265 1120
F: +27 11 265 1189

ADVANCED STERILIZATION PRODUCTS

Division of Medos International Sàrl

a **Johnson & Johnson** company