

ASP GLOSAIR™ 400 SYSTEM

Large-space decontamination.
Safe. Effective. Convenient.



ASP GLOSAIR™ 400 SYSTEM

Consistently protects your environment against pathogens.

The ASP GLOSAIR™ 400 SYSTEM delivers a mist of hydrogen peroxide and silver cations. Hydrogen Peroxide is an oxidizing agent, which kills microorganisms via the hydroxyl free radicals. Silver cations inactivate microorganisms by inhibition of protein synthesis.

- For patient wards, emergency rooms, operating rooms, and infectious disease units
- For room decontamination, just roll it in place
- Simple and quick; keeps your facility working

For use with ASP GLOSAIR™ 400 Cartridge.



TECHNICAL SPECS

| | |
|--------------------------------------|--|
| Cycle management | Track and trace capability |
| Mist particle size | 8 to 12 µm |
| Mean misting rate | 30 ml per minute ± 20% |
| Average decontamination ratio | 6ml/m ³ |
| Room size/ Volume treatment | 10 m ³ to 200 m ³ |
| Dimensions (mm) | 1080 H x 512 W x 620 D |
| Net weight | 48.8 kg |
| Power supply | 220-240 V / 50 Hz-7A (1680 W) 230-240 V / 60 Hz-7A (1680 W) |

COMPLIANCE DATA

| DIRECTIVES | HARMONIZED STANDARDS |
|--------------------|---|
| 2006/42/EC | NF EN 61010-1 |
| 2002/96/EC | Electrical and Electronic Equipment Directive |
| 2002/95/EC | ROHS Directive |
| 2006/66/EC | Batteries & Accumulators Directive |
| 2004/108/EC | EN 55022 CEI 61000-4-2 CEI 61000-4-3 CEI 61000-4-4 CEI 61000-4-5 CEI 61000-4-6 CEI 61000-4-11 CEI 61000-3-2 CEI 61000-3-3 |

Use biocides safely. Always read the label and product information before use.

FOR MORE INFORMATION

Contact for UK

ASP c/o J&J Medical Ltd.
Pinewood Campus,
Nine Mile Ride, Wokingham
Berkshire RG40 3EW,
England
T : +44 1 344 864 195
F : +44 1 344 871 171

Contact for Ireland

ASP c/o J&J Medical Ireland
Airtown Rd, Tallaght
Dublin 24, Ireland
T : +353 1 466 5200
F : +353 1 466 5340

Contact for Egypt

Johnson & Johnson
Medical Egypt
Florida Mall 5th Floor
1229 Square El Sheikh
Ali Gad
El Hak St. Heliopolis
Cairo, Egypt
T : +202 2268 5026
F : +202 2268 4674

Contact for Middle East

ASP c/o J&J Medical
Middle East
Dubai Healthcare City,
J&J Building, 3rd Floor
PO Box 505080, Dubai,
United Arab Emirates
T : +97 1 4 429 7235
F : +97 1 4 429 7250

Contact for South Africa

ASP c/o J&J Medical
Pty Ltd SA
PO Box 273, Midrand
Halfway House 1685
Republic of South Africa
T : +271 1 265 1120
F : +271 1 265 1189

Visit www.aspij.com/emea

ADVANCED STERILIZATION PRODUCTS

Division of Medos International Sàrl
a **Johnson & Johnson** company