



TECHNICAL SERVICES

STERRAD® 100NX® SYSTEM PLANNED MAINTENANCE ACTIVITIES PERFORMED DURING PM

DESCRIPTION	PM1	PM2
Drain the Vacuum Pump Oil	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Oil Mist Filter Element Replacement	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Air Filter/Fans	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Chamber Plastics	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
PM Counter	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Remove and Replace the Transfer and Transition Valves	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Converter Replacement	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
HEPA Vent Valve Filter	<input type="checkbox"/>	<input checked="" type="checkbox"/>
H ₂ O ₂ Monitor Optic Windows	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Hydrogen Peroxide Monitor Lamp	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Door Seal O-Ring Replacement	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Replace Thermistors	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Remove and Replace the Delivery Valve and Delivery Tubes	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Remove and Replace the Delivery Valve Air Filter	<input type="checkbox"/>	<input checked="" type="checkbox"/>



Follow us on Twittercom/ASPJJ



Like us on Facebook.com/ASPJJ



888.783.7723

ADVANCED STERILIZATION PRODUCTS

Division of Ethicon, Inc.

©Ethicon, Inc. 2013. All rights reserved.

a **Johnson & Johnson** company

AD-130044-01-US_B