



TECHNICAL SERVICES

STERRAD® 200 SYSTEM PLANNED MAINTENANCE ACTIVITIES PERFORMED DURING PM



DESCRIPTION	PM1	PM2
Change Vacuum Pump Oil	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Replace the Exhaust Filter	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Remove and Clean the Vaporizer Bowl	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Drain Water From the Air Storage Tank	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Drain Pressure Regulator Moisture Trap	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Clean Air Filters/Panel Assemblies	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Inspect/Clean Electrode (Replace after 2 years)	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Inspect/Clean Door, and Inside Chamber	<input checked="" type="checkbox"/>	<input checked="" type="checkbox"/>
Replace the Vaporizer Bowl O-Ring	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Baratron Zero Shift Test	<input type="checkbox"/>	<input checked="" type="checkbox"/>
H ₂ O ₂ Lamp Replacement	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Replace the Catalytic Converter	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Clean Platen Surface	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Replace Door Seal O-Ring	<input type="checkbox"/>	<input checked="" type="checkbox"/>
Replace Shelves, Stops, Rails, and Electrode Spacers	<input type="checkbox"/>	<input checked="" type="checkbox"/>



Follow us on Twittercom/ASPJJ



Like us on Facebook.com/ASPJJ



888.783.7723

ADVANCED STERILIZATION PRODUCTS

Division of Ethicon, Inc.

a **Johnson & Johnson** company