



THE ULTIMATE IN INSTRUMENT PROTECTION. Protecting your valuable, delicate surgical instruments and devices during sterilization, storage, and transport is vital. That's why Advanced Sterilization Products offers a full line of trays, silicone mats and instrument holders to secure and protect your instrument sets. APTIMAX® Instrument Trays feature a rugged, patented material that is compatible with most commonly used sterilization systems. It's easy to see why APTIMAX® Instrument Trays and Accessories offer the ultimate in protection for your valuable instruments.



FEATURES AND BENEFITS

Protection

- Durable material stands up to washer and sterilization processes

Standardization

- Fully validated for use in all STERRAD® Systems, pre-vacuum steam and ethylene oxide sterilization
- Full line of trays, mats and instrument holders allow you to customize instrument sets

Convenience

- Lightweight design and snug-fitting lid for ease in handling
- Stackable design to optimize space efficiency for storage

ADVANCED STERILIZATION PRODUCTS

Division of Ethicon, Inc.

a *Johnson & Johnson* company

APTIMAX® Instrument Trays & Accessories


APTIMAX® INSTRUMENT TRAYS

CODE	PRODUCT DESCRIPTION	DIMENSIONS	PACKAGING
13826	APTIMAX® Instrument Tray	12" x 4" x 2" (305 mm x 102 mm x 51 mm)	2/case
13827	APTIMAX® Instrument Tray	18" x 4" x 2" (457 mm x 102 mm x 51 mm)	2/case
13828	APTIMAX® Instrument Tray	24" x 4" x 2" (610 mm x 102 mm x 51 mm)	2/case
13829	APTIMAX® Instrument Tray	12" x 8" x 2" (305 mm x 203 mm x 51 mm)	2/case
13830	APTIMAX® Instrument Tray	18" x 8" x 2" (457 mm x 203 mm x 51 mm)	2/case
13831	APTIMAX® Instrument Tray	24" x 8" x 2" (610 mm x 203 mm x 51 mm)	2/case
13832	APTIMAX® Instrument Tray	7" x 3" x 1.2" (178 mm x 76 mm x 31 mm)	2/case
13833	APTIMAX® Instrument Tray	10" x 6" x 1.2" (254 mm x 152 mm x 31 mm)	2/case
13834	APTIMAX® Instrument Tray	15" x 6" x 1.2" (381 mm x 152 mm x 31 mm)	2/case
13835	APTIMAX® Instrument Tray	15" x 10" x 1.2" (381 mm x 254 mm x 31 mm)	2/case
13836	APTIMAX® Instrument Tray	11" x 10" x 4" (279 mm x 254 mm x 102 mm)	2/case
13837	APTIMAX® Instrument Tray	23" x 11" x 4" (584 mm x 279 mm x 102 mm)	2/case
13838*	APTIMAX® Instrument Tray	24" x 13" x 4" (610 mm x 330 mm x 102 mm)	1/case

* Not designed to fit in STERRAD® NX™ Systems



APTIMAX® Instrument Tray




APTIMAX® Instrument Trays & Mats

STERRAD® INSTRUMENT TRAY MATS

CODE	PRODUCT DESCRIPTION	DIMENSIONS	PACKAGING
99205	STERRAD® Instrument Tray Mat	6.5" x 2.5" (165 mm x 64 mm) fits tray 13832	2/case
99206	STERRAD® Instrument Tray Mat	9" x 5.5" (229 mm x 140 mm) fits tray 13833	2/case
99207	STERRAD® Instrument Tray Mat	14" x 5.5" (356 mm x 140 mm) fits tray 13834	2/case
99208	STERRAD® Instrument Tray Mat	14" x 9.5" (356 mm x 241 mm) fits tray 13835	2/case
99209	STERRAD® Instrument Tray Mat	11" x 7.5" (279 mm x 191 mm) fits tray 13829	2/case
99210	STERRAD® Instrument Tray Mat	17" x 7.5" (432 mm x 191 mm) fits tray 13830	2/case
99211	STERRAD® Instrument Tray Mat	23" x 7.5" (584 mm x 191 mm) fits tray 13831	2/case
99212	STERRAD® Instrument Tray Mat	10" x 10" (254 mm x 254 mm) fits tray 13836	2/case
99213	STERRAD® Instrument Tray Mat	21.5" x 10" (546 mm x 254 mm) fits tray 13837	2/case

APTIMAX® INSTRUMENT TRAY HOLDERS

CODE	PRODUCT DESCRIPTION	DIMENSIONS	PACKAGING
99400	APTIMAX® Instrument Tray Holder	0.4" x 0.8" x 0.4" (10 mm x 19 mm x 10 mm)	4/case
99401	APTIMAX® Instrument Tray Holder	0.8" x 0.4" x 0.8" (19 mm x 10 mm x 19 mm)	4/case
99402	APTIMAX® Instrument Tray Holder	0.2" x 0.8" x 0.2" (5 mm x 19 mm x 5 mm)	4/case
99403	APTIMAX® Instrument Tray Holder	0.8" x 0.2" x 0.8" (19 mm x 5 mm x 19 mm)	4/case
99404	APTIMAX® Instrument Tray Holder	0.4" x 0.4" x 0.4" (10 mm x 10 mm x 10 mm)	4/case
99405	APTIMAX® Instrument Tray Holder	0.6" x 0.6" x 0.6" (15 mm x 15 mm x 15 mm)	4/case
99406	APTIMAX® Instrument Tray Holder	0.2" x 0.2" x 0.2" (5 mm x 5 mm x 5 mm)	4/case
99407	APTIMAX® Holding Bar		4/case



APTIMAX® Instrument Tray Holder

FOR MORE INFORMATION

Contact an ASP representative at 1-888-STERRAD or visit www.aspjj.com.
For inquiries outside the U.S., please contact your local Advanced Sterilization Products representative.